

## Laboratory Genetic Metabolic Diseases

### Protocol for collecting CSF samples for metabolic screening

#### *Includes the following analyses:*

Neurotransmitter metabolites, 5-methyltetrahydrofolate (5-MTHF), amino acids, GABA, homocarnosine, B6 vitamers, pipecolic acid and pterins.

#### *General considerations*

- Preferably collect CSF samples between 8 and 10 h in the morning, before any medication is taken.
- Always mention the actual medication on the request form.
- Collect the CSF samples in small clean tubes without additives.
- Note the fraction numbers on the CSF samples. Reference values for neurotransmitter metabolites are determined for this specific fraction.

#### *Sample collection and storage*

Collect fractions in numbered tubes according to the following scheme:

Fraction	Volume	Instruction	Analysis
1	1 ml	On ice	Amino acids + special requests
2	1 ml	On ice	Neurotransmitter metabolites, 5-MTHF, pterins
3	0.5-1 ml	On ice	GABA, homocarnosine, pipecolic acid, B6 vitamers
4	0.5-1 ml	On ice and deproteinized	Lactate, pyruvate

- For newborns collect at least 2 fractions of 1 ml.
- Transport the CSF samples on ice to the lab, as fast as possible.
- In case of blood contamination immediately centrifuge samples prior to freezing and transfer the clear supernatants to labeled clean tubes.
- For lactate and pyruvate measurements deproteinize the sample according to protocol, see [www.amc.nl/lgmd/forms](http://www.amc.nl/lgmd/forms)
- **Freeze samples immediately at -80 °C and store until shipment** (if not possible, freeze at minimally -20 °C).

#### *Shipment*

- Send the frozen samples on sufficient dry ice to our laboratory by (international) courier together with a completed test request form Metabolite Diagnostics (available at [www.labgmd.nl](http://www.labgmd.nl)).